

SE-400

Electrosurgical Equipment



- High definition display
- Automatic detection function
- Memory function
- Detects monopolar and bipolar neutral plates
- Automatic compensation



SE-400

Functions:

- High definition digital display
- The machine has automatic detection function.
- It has a memory function, restart the emergence of the last used power setting machine
- It has a neutral electrode contact monitoring circuit, detection, anad patient contact plate area can be effective contact area, evaluation of plates and skin system testing, once found the contact area is reduced to a dangerous level, automatically cut of the output system and alarm
- It can detect monopolar and bipolar neutral plates
- It has automatic compensation function, stable output power, makes the loss to a minimum
- It can be used for under water cutting, including lipid rich tissue cutting area
- It can be use with the endoscope and hysteroscopy

Specifications:

Output power:	400W
Five output modes:	Cutting (400 W) Blend Cut (150 W) Soft Coag (100 W) Strong Coag (80 W) Bipolar Coag (50 W)
Atmospheric pressure range:	86.0 ~ 106.0 kPa;
Power supply:	250V±22V, 50Hz±1
Working frequency:	360 ~ 460 kHz;
The rated power of the equipment:	880VA±10%;
Type of device:	CF;
Machine power consumption:	≤1000VA;

SANTAMED

ZHENGZHOU SANTA MEDICAL CO.,LTD

Room 301 Building-6, Jincheng Time Square, Jinshui Zone, Zhengzhou450000 Henan Province PR,China

+ 86-177 19820906

sales@santamedhk.com; cloverye@santamedhk.com